



Medical Device Details

Sponsor: CARSL Consulting
P O Box 766
Hastings
New Zealand

WAND Reference: 161208-WAND-6N7UFT

Sponsors Own Reference: XCITE CLINICAL STATION

GMDN: Stimulator, neuromuscular - specify [35370]

Class: IIa

Intended Purpose: Xcite is intended to be used for: 1. relaxation of muscle spasms; 2. prevention or retardation of disuse atrophy; 3. increasing local blood circulation; 4. muscle re-education; 5. immediate post-surgical stimulation of calf muscles to prevent venous thrombosis; 6. maintaining or increasing range of motion.

Manufacturer: Restorative Therapies Inc 8098 Sandpiper Circle, Suite M
Nottingham
Maryland
United States

ARTG ID:
The device is supplied unsterilised.
The device is not intended to be invasive.
The device is not intended for single use.
The device is an active device.
The device does not contain material or ingredients of microbial origin.
The device does not contain material or ingredients of recombinant origin.
The device does not contain material or ingredients manufactured or formulated using a genetically modified organism.
The device does not contain material or ingredients of human origin.
The device does not contain human blood or its components.

The device consists of: System
 Components:
Component: Cart, optional cart
Component: Database, RTILink.com online database
Component: Electrodes, Adhesive electrodes
Component: Xcite Stimulator, Xcite Stimulator with: - one or two 6 channel stimulators - an easy-to-use touch interface - backup battery power - AC/DC power a
The device does not contain material or ingredients of animal origin.
The device is not medicated.
The device is not formulated.
The product does not contain a medicine that has consent for marketing in New Zealand.
The product does not contain a medical device which incorporates a medicine as an integral part and that has an action ancillary to the device.

Status: Active